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
Three-points Meter

Activity, food intake and water intake monitor

<amount of activity, food intake and water intake / biological rhythm >

※Multi chamber system (1 – 16 chamber)



Special cages designed to house up to four subjects are utilized in this test. Each cage is equipped with a CCD camera placed directly above and captures images at 1 - 2 frames/sec. To reduce any interruptions during image capture, both water bottle and food box are fixed on the sides of the cage. When only analyzing a single subject within the cage, each image frame is converted to a binary image and the elliptical subject particle is tracked by using the center of the particle. With this analysis, the software can measure the subject's distance traveled every minute. When two to four subjects are housed in the cage, each image is converted to a binary image and the software calculates the difference of changed area from the previous image to measure the amount of activity. Furthermore, to analyze whether the subjects come in contact with another during testing, the software analyzes the number of particles in each image frame. This analysis shows when the subjects are in contact with others or are separated. If a large capacity hard disk is connected to the computer, the software can also store serial images during testing. The software automatically analyzes and stores all data files and images throughout testing.

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Activity data acquisition software

Food and Water intake data acquisition software

Food intake and water intake are calculated by the difference between the measured weights. The apparatus is designed to eliminate dropped food, bedding material, and defecations to be added into the measurements. By using the software, time series data can automatically be saved as text files.

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This system can measure for 3 kinds of amounts: 1 ambulatory activity, 2 food intake and 3 water intake, at the same time within one cage. It is possible to measure at the same time such amount within 16 cages at the maximum (16 individuals) by means of a computer. For measuring activity, the position of the target is detected by interrupting the infrared beams positioned in the X and Y axis around the cage. This device detects the movement in position of the whole body of the target animal, without being misled by the movement of the tail and any other body. The position and movement of the target within the cage are automatically measured. The food box and the water bottle are placed under a precision scale and the software periodically measures the weight of the box.